

## Phosphorylated Compounds Analyzed with Cogent Diamond Hydride Columns - Tips and Suggestions

Date: 14-JANUARY-2016 Last Updated: 30-JUNE-2025

Analyzing phosphorylated compounds by HPLC can be challenging due to their high polarity and tendency to interact with metal surfaces. However, the Cogent Diamond Hydride™ column, used in Aqueous Normal Phase (ANP) mode, offers a powerful solution for retaining and separating a wide range of phosphate-containing analytes.

### **Compatible Compound Types**

The column has demonstrated excellent retention and separation for:

- Phosphorylated sugars
- Herbicides
- Nucleotides
- Isoprenoid phosphates and other phosphate-functionalized molecules

### **LC-MS Detection Tips**

When using LC-MS, monitor for  $[M - H]^-$  ions in negative ion mode extracted ion chromatograms (EICs) for optimal detection of phosphorylated species.

### **Recommended Mobile Phase Conditions**

Use the following mobile phase system for ANP:

- Solvent A: DI water + 10 mM ammonium acetate (or ammonium formate)
- Solvent B: 95:5 Acetonitrile:DI water + 10 mM ammonium acetate

Start with a high %B and apply a gradient to a moderate or low %B to adjust retention. Fine-tune the gradient based on your compound's behavior.

### **Troubleshooting Peak Shape Issues**

To improve peak shape and reproducibility, consider these strategies:

1. Adjust pH: Add a small amount of dilute ammonia to raise the mobile phase pH to ~7.0.
2. Avoid Sodium Contamination: Use Teflon bottles instead of glass to prevent sodium leaching.
3. Metal Ion Interference: Add 5–10  $\mu$ M EDTA to the mobile phase to chelate trace metals that can distort peaks of anionic compounds.

4. Column Hardware: Use metal-free, coated HPLC columns to avoid unwanted interactions with stainless steel surfaces.

## Ordering Information

Click [HERE](#) for ordering the Cogent Diamond Hydride™ column or to learn more about its specifications and available formats, please visit our [product page] or contact our technical support team.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)