

## How does ammonium fluoride affect chromatography using Cogent TYPE-C HPLC columns - FAQ

DATE: 24-OCTOBER-2016 Last Updated: 6-OCTOBER-2025

Ammonium fluoride offers several advantages when used as a mobile phase additive in LC-MS applications when used with Cogent TYPE-C HPLC Columns. The columns have a **longer lifetime** than standard HPLC columns due to the silica hydride particle that is somewhat resistant to the corrosive effects of the acid.

- **Acts as a buffer**, helping to control pH
- **Enhances peak shapes**, leading to better chromatographic performance
- **Improves MS sensitivity** for certain analytes

For more detailed information, refer to the following third-party journal article: J. Pesek, M. Matyska, "Ammonium fluoride as a mobile phase additive in aqueous normal phase chromatography," J. Chromatogr. A. 1401, 2015, pp. 69–74.

### Important Note:

*Ammonium fluoride is corrosive to glass, so Teflon containers are recommended. Additionally, concentrations above 1 mM may cause irreversible damage to HPLC columns. While Cogent TYPE-C™ Silica columns have been reported to exhibit greater durability compared to conventional columns, ammonium fluoride may still reduce their lifespan to some extent.*



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)