

Selectivity Modification in an ANP Method Using Different Mobile Phase pH - Tips & Suggestions

Date: 15-OCTOBER-2012 Last Updated: 30-NOVEMBER-2025

Optimizing Selectivity in Aqueous Normal Phase (ANP) HPLC

One of the most effective and straightforward ways to modify selectivity in ANP is by adjusting the mobile phase pH.

Why pH Matters

The ionization state of both the analytes and the stationary phase plays a critical role in ANP retention. Changing the pH can significantly impact chromatographic behavior.

Key Considerations

- Use mobile phase additives that are compatible with your stationary phase and system.
 - Avoid phosphates in LC-MS applications.
- Common additives include:
 - 0.1% formic acid or acetic acid
 - 10 mM ammonium acetate or ammonium formate
- Refer to the Cogent™ column specifications for allowable pH ranges, especially for Diamond Hydride™ columns.

Solvent Choice

The organic component in ANP methods also influences retention. For example:

- Acetone vs. Acetonitrile can produce markedly different retention profiles.
- Combining both solvents may help optimize selectivity:
- A peak pair that co-elutes with acetone may separate with acetonitrile—and vice versa.
- A mixed solvent approach can sometimes achieve complete resolution of all peaks.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com