

## Modify Selectivity in an ANP Method Using Different Mobile Phase pH - Tips & Suggestions

*Date: 15-OCTOBER-2012 Last Updated: 30-NOVEMBER-2025*

### Optimizing Selectivity in Aqueous Normal Phase (ANP) HPLC

One of the most effective and straightforward ways to modify selectivity in ANP is by adjusting the mobile phase pH.

#### Why pH Matters

The ionization state of both the analytes and the stationary phase plays a critical role in ANP retention. Changing the pH can significantly impact chromatographic behavior.

#### Key Considerations

- Use mobile phase additives that are compatible with your stationary phase and system.
  - Avoid phosphates in LC-MS applications.
- Common additives include:
  - 0.1% formic acid or acetic acid
  - 10 mM ammonium acetate or ammonium formate
- Refer to the Cogent™ column specifications for allowable pH ranges, especially for Diamond Hydride™ columns.

#### Solvent Choice

The organic component in ANP methods also influences retention. For example:

- Acetone vs. Acetonitrile can produce markedly different retention profiles.
- Combining both solvents may help optimize selectivity:
- A peak pair that co-elutes with acetone may separate with acetonitrile—and vice versa.
- A mixed solvent approach can sometimes achieve complete resolution of all peaks.

Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)