

Can I backflush the Cogent column filter for longer in-service life - FAQ

Date: 3-DECEMBER-2022 Last Updated: 11-OCTOBER-2025

Can I Backflush the Cogent Column Filter to Extend Its Service Life?

While the Cogent column filter is technically capable of being backflushed, ***we do not recommend this practice.***

Backflushing may cause the internal filter element to shift or become dislodged, even if no visible clog is present. A misaligned filter can compromise its effectiveness, leading to inconsistent flow or reduced protection for your column.

For this reason, backflushing is not supported as a maintenance method for Cogent column filters.

👉 Click [HERE](#) for ordering information and to view images of the Cogent column filters.

