

## Cogent Column Filter Backflushing is Not Recommended - Tech Information

*Date: 3-DECEMBER-2022 Last Updated: 11-OCTOBER-2025*

**While the Cogent column filter is technically capable of being backflushed, we *do not* recommend this practice .**

Backflushing may cause the internal filter element to shift or become dislodged, even if no visible clog is present. A misaligned filter can compromise its effectiveness, leading to inconsistent flow or reduced protection for your column.

For this reason, backflushing is not supported as a maintenance method for Cogent column filters.

👉 Click [HERE](#) for ordering information and to view images of the Cogent column filters.



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)