



## Thiopurines Analysis by LC-MS - AppNote

### Separation using LCMS Gradient Method

The Cogent Diamond Hydride Column with gradient elution was used for analysis of Thiopurines. *Figure A* shows the Peak of 6-TG, which also can be analyzed using Isocratic conditions. Good retention and symmetrical peak shape were obtained under the analysis conditions.

*Figure B* represents two separated Thiopurines (Mono and Di- Phosphate forms). *Figure C* shows two Inosine compounds, one with an additional Methyl group, being separated.





### Peaks:

A: Thioguanine (6-TG) at  $m/z = 168.0338 [M+H]^+$

B: 6-Thioguanosine -5'-Phosphate (TGMP) at  $m/z = 380.3$ , 6-Thioguanosine -5'-Diphosphate (TGDP) at  $m/z = 460.3$

C: 6-Methyl-Thioinosine-5'-Monophosphate (meTIMP) at  $m/z = 379.3$  and 6-Thioinosine-5'-Triphosphate (TITP) at  $m/z = 525$

### Method Conditions

**Column:** Cogent Diamond Hydride™, 4  $\mu\text{m}$ , 100  $\text{\AA}$

**Catalog No.:** [70000-15P-2](#)

**Dimensions:** 2.1 x 150 mm

### Mobile Phase:

- A: DI Water / 50% Methanol / 0.1% Formic Acid (v/v)
- B: Acetonitrile / 0.1% Formic Acid (v/v)

### Gradient:

Time (minutes)	%B
0	100
12	30
14	30
15	0
19	0
20	100

**Post Time:** 2 minutes

**Flow rate:** 0.4 mL / minute

**Detection:** ESI – POS - Agilent 6210 MSD TOF Mass Spectrometer

**Injection vol.:** 1  $\mu\text{L}$

**Sample Preparation:** 0.4 mg / mL solutions in DI Water. For MS analysis, samples were diluted 1:100 into 50% Acetonitrile / 50% DI Water mixture. Before injection, samples were filtered through a 0.45  $\mu\text{m}$  Nylon Syringe Filter (MICROSOLV Tech Corp.).



**Attachment No 322 Analysis of Thiopurines pdf 0.3 Mb** [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)