

Explanation of Wavelength and Absorbance Precision Test of the UV-vis Software -Tech Information

DATE: 27-JUNE-2025 Last Updated: 28-JUNE-2025

Wavelength Calibration

When working with spectrophotometers that use mechanical monochromators (devices that select specific wavelengths of light), it's important to check how precisely the instrument can return to a chosen wavelength—this is called wavelength precision or reset capability.

Here's how the Wavelength Precision Tests (resetability) test is done in our UV-vis Software:

1. Set the wavelength: The wavelength resetability can be performed at any nominal wavelength that the user chooses. In the example, it is set at 278 nm.
2. Read the absorbance value.
- 3.
4. Change the wavelength setting to several nm above or below the nominal wavelength, then return to the nominal wavelength.
5. Read the absorbance value.
6. Repeat this six times (do three readings above, and three below the nominal wavelength, alternating).
7. The Standard Deviation (SD) and Coefficient of Variation (CV) are calculated from those absorbance values—this tells you how consistent and precise the instrument is at resetting the same wavelength.

<u>Wavelength and Absorbance Precision Tests:</u>			
Wavelength for Precision (Resetability) Test:		278 nm	
Trial No.	Absorbance @ Wavelength	λ Precision Results:	
1		Average:	#DIV/0!
2		Min:	0.000
3		Max:	0.000
4		Std. Deviation:	#DIV/0!
5		%RSD (CV):	#DIV/0!
6			

Absorbance precision, as shown in the right-hand table, involves removing and reinserting the cuvette six times, with absorbance measured after each insertion.

Absorbance Solution for Precision Test:

Solution No. 2

Trial No.	Absorbance Solution:	Photometric Precision	
		Results:	
1		Average:	#DIV/0!
2		Min:	0.000
3		Max:	0.000
4		Std. Deviation:	#DIV/0!
5		%RSD (CV):	#DIV/0!
6			

It's important to note that the two tables represent separate measurements and are not directly related.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.comWebsite: www.mtc-usa.com