



Ribose and Xylose Analyzed by HPLC ELSD - App Note

Date: 10-OCTOBER-2025 Last Updated: 10-OCTOBER-2025

Sugars are challenging to analyze by HPLC due to their polarity and lack of UV absorbance. ELSD enables detection without relying on chromophores. While amine columns are often used for sugar retention, they can react with aldehydes, reducing column life by formation of Schiff bases, resulting in irreversible deactivation of the ligand's retention functionality. The Cogent Amide Column avoids this issue with a less reactive ligand, offering stable retention and separation of Ribose and Xylose.

See this method below for great separation of these two epimers.



Peaks:

1. D-Ribose
2. D-Xylose

Method Conditions:

Column: Cogent Amide™, 4 µm, 100 Å

Catalog No.: 40036-10D

Dimensions: 4.6 x 100 mm

Mobile Phase: 95% Acetonitrile / 5% DI Water / 0.1% Triethylamine (TEA) (v/v)

Flow Rate: 1.0 mL / minute

Detection: ELSD (Gain 8, Temp 50C, Nitrogen: 3.5lb)

Injection Volume: 3 µl

Sample Preparation: D-Ribose and D-Xylose reference standards (1 mg/mL) in diluent of 50% Acetonitrile / 50% DI Water / 0.1% TEA (v/v)

Note: Ribose and Xylose are aldopentoses that differ only by a chiral center. In addition to the open chain forms, these sugars exist in equilibrium with ring forms (five or six membered) as well as α and β anomers. Both sugars are highly polar and not generally suitable for conventional Reversed Phase retention.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com