


Creatine and Creatinine Analyzed with ELSD and UV - AppNote

Date: 13-OCTOBER-2025 Last Updated: 13-OCTOBER-2025

Creatine and creatinine are key biomarkers in clinical diagnostics. Their accurate separation is essential for reliable analysis. This application note describes an HPLC method using ELSD and UV detection optimized for the efficient separation and quantification of creatine and creatinine.



Peaks:

- 1. Creatinine**
- 2. Creatine**

Method Conditions:

Column: Cogent Diamond Hydride™, 4 µm, 100 Å

Catalog No.: 70000-7.5P

Dimensions: 4.6 x 75 mm

Mobile Phase: 70% Acetonitrile / 30% DI Water

Flow Rate: 1.0 mL / minute

Detection:

A: UV 210 nm

B: ELSD (Gain 8, Temp 50°C, Nitrogen: 3.5 lb)

Injection Volume: 2 µL

Sample Preparation: Creatinine (0.45 mg/mL) and Creatine (0.26 mg/mL) reference standards in diluent of 50% Acetonitrile / 50% DI Water

Note: Creatine is stored in muscles and used for quick energy during intense activity, while creatinine is a waste product formed when creatine breaks down—making it a key marker for kidney function in medical tests.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com