

Are Polypropylene Vials Suited to Avoiding Adsorption Issues with Proteins Compared to Glass - FAQ

Date: 10-JUNE-2013 Last Updated: 27-OCTOBER-2025

Yes, polypropylene as a vial material is a good choice for protein analysis because normal borosilicate glass vials are known to adsorb some types of analytes such as proteins, peptides, and basic compounds such as amines.

This leads to lower ~~recovery~~ and is caused by silanols on the glass surface. Polypropylene does not have these groups and therefore adsorption is not an issue with these vials.

It is important to note that some plastic vials are made in high speed machines and contain a lot of "mold release agents". These vials would not be suitable for many proteins as the extractables could cause damage to your samples. Also, these agents could also cause erroneous peaks in LCMS. The Mass Spec compatible vials are most suitable as they are not made in high speed machines and we avoid the use of "mold release agents" making them more suitable for proteins and minimizing extractables.

Click [HERE](#) for ordering information and pictures of LCMS compatible plastic vials.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com