

## Considerations for keeping gases dissolved in solution when using sparging stones - Tips & Suggestions

**Question:** I am using the sparging stones (*MTC catalog# 48222-06*) in my production lines feeding nitrogen into solution. After switching to a bigger dish with more volume, I noticed that the concentration of nitrogen does not reach the same level anymore. What can you suggest for a sparging stone that can handle higher flow rates (*current method around 40mL/min*) ? With our current settings, the nitrogen just bubbles to the surface.

**Answer:** Here are some potential solutions to the issue:

1. Is the stone mounted in the dish?
2. Can an agitator be added? This may make the process of getting nitrogen into the solution more effective as 40 mL/min is a very low flow rate.
3. Can a mesh be added over the stone? This may help keep the gas in solution.

Click [HERE](#) for ordering information and pictures of sparging stones.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)