



## Pyrilamine and 4-Amino-3-Chloropyridine Analysis with HPLC - AppNote

### Unique Selectivity on a Cogent Amide Stationary Phase

The Cogent Amide column offers unique selectivity that may not be readily attainable with other phases. Two test solutes shown in this application note (*Pyrilamine* and *4-Amino-3-Chloropyridine*) were baseline separated on the Cogent Amide column (*Figure A*), but they co-eluted with no resolution on a different Cogent column using otherwise equivalent method conditions (*Cogent Diamond Hydride™*, *Figure B*). The presence of the Amide ligand provides additional selectivity that can make a significant difference in resolving closely-eluting compounds such as these.



#### Peaks:

1. Pyrilamine

## 2. 4-Amino-3-Chloropyridine

### Method Conditions

**Column:** Cogent Amide™, 4µm, 100Å

**Catalog No.:** [40036-05P](#)

**Dimensions:** 4.6 x 50mm

### Mobile Phase:

A: 90% DI Water / 10% Acetonitrile / 0.1% Formic Acid (v/v)

B: B: Acetonitrile / 0.1% Formic Acid (v/v)

### Gradient:

| Time<br>(Minutes) | %B |
|-------------------|----|
| 0                 | 90 |
| 1                 | 90 |
| 7                 | 50 |
| 8                 | 90 |

**Post Time:** 3 minutes

**Flow rate:** 1.0 mL/minute

**Detection:** UV 244 nm

**Injection vol.:** 2µL

### Sample Preparation:

100 mg/L Pyrilamine and 4-Amino-3-Chloropyridine reference standards in diluent of 50/50 solvent A/solvent B. Peak identities confirmed with individual standards.

---

**Note:** Amine-containing compounds such as Pyrilamine and 4-Amino-3-Chloropyridine can be difficult to analyze using conventional silica- based stationary phases. These columns have residual silanol groups on the surface that can interact electrostatically with Amines, causing peak tailing. Chromatographers use a number of strategies to avoid these issues, such as use of ion pair agents or endcapping. However, Cogent TYPE-C Silica phases are virtually free of silanols, and therefore good peak shapes can be obtained without these workaround method strategies.

---



**Attachment No 352 Pyrilamine and 4-Amino-3-Chloropyridine.pdf** 0.4 Mb [Download File](#)