

## Do You Really Need High pH for Retention in HPLC - Tips & Suggestions

*Date: 25-FEBRUARY-2013 Last Updated: 30-NOVEMBER-2025*

Traditionally, some analyses have required very high mobile phase pH (up to 10–12) to achieve acceptable retention for compounds with high pK<sub>a</sub> values. However, using extreme pH is not always necessary—and it comes with significant drawbacks.

### Risks of High pH

- **System Damage:** Prolonged exposure to high pH can degrade pump seals, tubing, fittings, and detector cells.
- **Column Issues:** At very high pH and moderate ionic strength, silica dissolution can occur, shortening column life.

---

### Alternative Approach

With Cogent™ TYPE-C™ HPLC Columns and Aqueous Normal Phase (ANP), many compounds can be retained effectively within the typical HPLC pH range of 2–7.

- Some methods still require high pH for analyte solubility—if so, maintain high pH.
- In ANP, the goal is to fully ionize compounds, not neutralize them.

---

### Rule of Thumb

- Acidic analytes → Basic mobile phase
- **Basic analytes → Acidic mobile phase**

