

Difference Between Waters-Type Fittings and Standard 10-32 Fittings in HPLC - HPLC Primer

Date: 2-MARCH-2024 Last Updated: 19-FEBRUARY-2026

Overview

Waters-type fittings and standard 10-32 fittings are both widely used in HPLC systems, yet they differ in design, seat geometry, and compatibility with instrument ports.

These differences impact sealing performance, tubing engagement, and whether a given fitting style can be interchanged safely. Understanding how each fitting works ensures leak-free, low-dead-volume chromatographic connections.

What Are Waters-Type Fittings?

Design Characteristics

- Waters-type fittings use a unique concave seat geometry, engineered specifically for Waters brand HPLC instruments and components.
- The internal geometry is designed to mate with Waters-style column end fittings, often incorporating a coned or chamfered sealing surface.
- They may require specific ferrule shapes that match the Waters port design.

Performance Notes

- Waters fittings ensure:
 - Proper tube centering
 - Optimal sealing against Waters-designed port geometry
 - Minimized dead volume in matching Waters hardware

Compatibility

- Waters-type fittings cannot always be substituted with generic 10-32 fittings because the seat angles and sealing surfaces differ.
 - Using a mismatched fitting can result in:
 - Leaks
 - Poor peak shape
 - Ferrule deformation
 - Damage to the receiving port
-

What Are Standard 10-32 Fittings?

Design Characteristics

- Standard 10-32 fittings use a universal 10-32 UNF thread and a coned (typically 60°) seat common across many HPLC systems.
- Compatible with a broad range of:
 - Columns
 - Unions
 - Adapters
 - Valves and detectors
 - Instrument brands that follow standard coned-port geometry

Performance Notes

- Provide excellent sealing when paired with correctly matched ferrules and coned ports.
 - Offer greater interchangeability than Waters-type fittings.
-

Key Differences Between Waters-Type & Standard 10-32 Fittings

Threading

- Both use 10-32 threads, but the seat geometry is different.

Port Geometry

- Waters: Concave or specialized seat design.
- Standard 10-32: Conventional coned seat.

Ferrule Compatibility

- Waters fittings require Waters-compatible ferrules.
- Standard 10-32 fittings use generic ferrules sized for 1/16" OD tubing.

Interchangeability

- Waters fittings should not be used in standard 10-32 coned ports.
- Standard 10-32 fittings should not be used in Waters-type ports.

Result of Mismatch

- Poor sealing
 - Increased dead volume
 - Leaks
 - Damage to the column's receiving port
 - Misalignment of tubing inside the fluidic channel
-

When to Use Each Fitting Type

Use Waters-Type Fittings When

- Connecting to Waters columns or components designed for this seat geometry
- Ensuring the tightest seal with Waters-engineered hardware
- Replacing original fittings on Waters instruments

Use Standard 10-32 Fittings When

- Working with instruments or columns from Agilent, Shimadzu, Thermo, and other brands using standard coned-port designs
- Building general HPLC systems with high interchangeability
- Connecting tubing to unions, tees, and adapters with standard 10-32 ports

Click [HERE](#) for information about these adapters that allow you to use standard fittings with Waters Columns or instruments.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com