

## UV Light Pass Does Not Pass Through PTFE Tubing Easily - Tech Information

Date: 10-October-2025 Last Updated: 20-FEBRUARY -2026

### Overview

PTFE (Polytetrafluoroethylene) is well-known for its exceptional resistance to chemical attack, thermal stress, and environmental exposure—including ultraviolet light. Its semi-crystalline structure allows it to reflect a significant portion of UV energy while transmitting some wavelengths without rapid degradation.

However, not all UV exposure is benign; the wavelength and intensity of the UV source determine whether PTFE remains stable or begins to degrade.

---

### UV Transmission and Wavelength Effects

#### 1. Transparency Across Most UV Wavelengths

PTFE shows general transparency across much of the UV spectrum, giving it a strong reputation for UV resistance.

#### 2. Strong Absorption Below 240 nm

At wavelengths below 240 nm, PTFE begins to absorb UV intensely, leading to:

- Surface chemical bond breakage
- Localized heating
- Photothermal degradation

These high-energy photons can alter the tubing surface, making it unsuitable for applications involving deep-UV or excimer-based illumination.

#### 3. No Degradation Above 400 nm

Visible and near-IR region exposure does not degrade PTFE, making it stable under typical laboratory lighting and higher wavelengths.

---

### Influence of Wall Thickness - Thinner Walls Allow More UV Transmission

Thin-walled PTFE tubing will transmit more UV—particularly at longer wavelengths—than thick-walled tubing.

However, even when transmitting UV, thin walls may still suffer localized damage if the UV energy is concentrated or intense.

---

## Degradation Under Intense UV Sources

PTFE's semi-crystalline microstructure can trap and scatter photons, making it vulnerable to:

- Localized photothermal degradation
- Surface morphology changes
- Structural weakening around irradiated zones

Excimer lasers and other high-intensity UV systems significantly increase degradation risk.

---

## Practical Considerations for Laboratory Use

PTFE tubing can be safely used in:

- General UV-exposed environments
- Fluorescence-based systems (non-deep-UV)
- Outdoor or ozone-producing areas

Avoid using PTFE tubing in:

- Deep-UV applications (<240 nm)
- High-intensity UV irradiation paths
- Environments with focused laser exposure

Click [HERE](#) for PTFE Tubing Ordering Information.

---

**NOTE:** While PTFE tubing is generally considered to be UV-resistant, exposure to intense UV sources such as excimer lasers can lead to localized photothermal degradation. This alters the surface morphology and is influenced by PTFE's semi-crystalline microstructure, which can trap photons and scatter light, affecting the extent and nature of the damage.

---

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)