

## Autoclavability of Stainless Steel Tubing and Fittings for HPLC Use - HPLC Primer

*Date: 12-MAY-2012 Last Updated: 15-FEBRUARY-2026*

### Introduction

Sterilization of HPLC components is sometimes necessary in bioanalytical, microbiological, or sterile-processing workflows. Stainless steel tubing and fittings—commonly made from 316 stainless steel—are routinely used in chromatography systems and must withstand sterilization procedures without compromising mechanical integrity or chromatographic performance.

Understanding which sterilization methods are compatible is essential for maintaining system reliability and preventing equipment damage.

---

### Compatibility of 316 Stainless Steel With Autoclaving and Other Sterilization Methods

316 stainless steel is widely recognized for its corrosion resistance, high-temperature stability, and chemical durability.

- Autoclaving: Many scientists have successfully sterilized 316 stainless-steel tubing and fittings using steam autoclave cycles.
  - Gamma irradiation: This non-thermal method is compatible with stainless-steel components, making it useful for sensitive systems requiring low-residue sterilization.
  - Ethylene oxide (EtO) sterilization: Also compatible, offering a cold-sterilization alternative for heat-sensitive assemblies while maintaining material integrity.
- These methods do not compromise the structural performance of 316 stainless-steel fittings or tubing when used appropriately.

---

### Practical Considerations for Sterilizing Stainless Steel Components

Although 316 stainless steel is highly resilient, proper handling ensures optimal results.

- Remove sensitive polymeric components: Ferrules, seals, or polymeric housings may not withstand high-temperature or radiation-based sterilization.
- Ensure complete drying after autoclaving: Residual moisture can promote corrosion in downstream system parts if not removed.
- Follow system manufacturer guidelines: Some HPLC instruments have specific recommendations for sterilization cycles or pressure limits.
- Avoid harsh chlorinated cleaners: While 316 stainless steel is corrosion-resistant, chloride exposure at elevated temperatures may accelerate pitting.

These precautions help maintain the longevity and performance of stainless-steel HPLC tubing and fittings.

---

## When Sterilization Is Appropriate in HPLC Workflows

Sterilization is typically needed when working with:

- Microbial or cell-based samples
- Sterile pharmaceutical preparations
- Biological matrices prone to contamination

- Environments requiring aseptic handling between runs

In these scenarios, ensuring the sterility of flow paths—including tubing and fittings—prevents biofilm formation, contamination artifacts, and instrument downtime.

---

## Conclusion

316 stainless-steel tubing and fittings used in HPLC systems are compatible with several sterilization methods, including autoclaving, gamma irradiation, and ethylene oxide treatment.

Their durability and corrosion resistance make them suitable for sterile analytical environments, provided that proper handling protocols are followed to protect system integrity and performance.

Click [HERE](#) for stainless steel tubing ordering information

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)