

Differences Between Cogent Amide Columns and Other Brands of Amide - Tech Information

Date: 7-JULY-2018 Last Updated: 16-FEBRUARY-2026

Overview

When choosing an amide-type stationary phase, subtle chemistry details can make or break your method—especially for carbohydrates and basic pharmaceuticals. The Cogent Amide™ column is engineered on TYPE-C™ silica and leverages direct silicon–carbon (Si–C) bonds.

This backbone enables robust operation in Reversed-Phase (RP), Normal Phase (NP), and Aqueous Normal Phase (ANP) environments with improved resistance to hydrolysis compared with many conventional amide columns that rely on hydrolyzable linkages.

Key Performance Advantages of Cogent Amide™

1. Wider Method Flexibility (RP / NP / ANP)

Cogent Amide performs across multiple chromatographic modes, allowing you to tune selectivity for amines, simple sugars (e.g., glucose, galactose, sucrose), and basic drugs (e.g., fluoxetine / Prozac®, tizanidine) without changing hardware.

2. Hydrolysis Resistance

The Si–C anchoring in TYPE-C silica significantly reduces the hydrolysis that often plagues traditional amide phases, supporting longer column lifetimes and more stable retention in aqueous or mixed mobile phases.

3. No Schiff Base Pitfalls

Cogent Amide does not form Schiff bases with aldehydes and is not deactivated by them—avoiding the semi-permanent column modification that can derail carbohydrate analyses.

4. Cleaner Carbohydrate Separations vs. “Amino” Columns

“Amino” columns (bearing basic, ionizable amine groups) can interact strongly and sometimes irreversibly with carbonyl compounds, altering surface chemistry over time. In contrast, amide functional groups are **non-basic and non-ionizable**, minimizing these complications and preserving reproducibility.

Amide vs. “Amino” Columns — Practical Implications

• Surface Chemistry

- *Amino columns*: Basic/ionizable; prone to strong interactions with aldehydes and formation of Schiff bases.
- *Amide columns (Cogent Amide)*: Neutral, non-ionizable; do not create semi-permanent aldehyde associations.

- Method Robustness
 - *Amino columns*: Risk of gradual, hard-to-reverse performance shifts when exposed to carbonyl compounds (e.g., reducing sugars).
 - *Cogent Amide*: Maintains consistent performance with sugars and amine-containing analytes, reducing method drift and revalidation cycles.
 - Analyte Scope
 - *Amino columns*: Often used historically for sugars but susceptible to deactivation artifacts.
 - *Cogent Amide*: Designed for sugars and basic pharmaceuticals with dependable retention across RP/NP/ANP.
-

Application Examples (Representative)

- Simple Carbohydrates: Glucose, galactose, sucrose — separated without column modification concerns.
 - Basic Pharmaceuticals: Fluoxetine (Prozac®), tizanidine — robust retention with mixed-mode flexibility.
-

Why TYPE-C™ Silica Matters

TYPE-C silica's direct Si–C linkage is the chemical foundation behind Cogent Amide's durability and versatility. It supports mixed-mode selectivity and resists hydrolysis, enabling the column to function predictably in aqueous, organic-rich, or gradient conditions common to RP, NP, and ANP workflows.

Takeaway

If your methods involve sugars, amines, or basic drugs, and you need mode flexibility without sacrificing stability, Cogent Amide™ provides a reliable alternative to traditional amide and especially to “amino” columns, minimizing Schiff base formation and hydrolysis-driven drift.

FAQs

Q1. Can I use Cogent Amide™ for both sugars and basic amines in one lab?

Yes. The phase is effective for both compound classes and operates in RP, NP, and ANP styles, simplifying platform standardization.

Q2. Why avoid “amino” columns for reducing sugars?

Because amino phases are basic and can form Schiff bases with aldehydes, leading to semi-permanent changes and deteriorating performance.

Q3. What's the advantage of TYPE-C™ silica here?

Si–C bonds confer hydrolysis resistance and enable mixed-mode retention control, supporting broader method design with better durability.

[Click HERE for Cogent Amide HPLC Column Ordering Information.](#)



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com