

Total Method Development Kit

Best HPLC Tool for Hydrophilic Compounds



- Save Money & Time
- 3 Modes Capable
- Wide pH and Temperature Range



A suite of columns for most mixtures of analytes will make finding and analyzing your hydrophilic and hydrophobic compounds easier and faster. Useful with HPLC or LCMS, the suite of columns gives you a range of retention power ranging from the most polar to complex mixtures of hydrophobic and hydrophilic compounds.

Add to your current favorite columns in your scouting, problem solving or routine strategies. Save time and get great reproducible results.

Each suite contains 1 each of the following HPLC Columns based on TYPE-C Silica™ Cogent Silica-C™, Cogent Diamond Hydride™, Cogent Bidentate C8™, Cogent UDC-Cholesterol™, Cogent Bidentate C18.

For more information:
www.mtc-usa.com/wcode_cog_suite.asp

Cat. No.	Description	
5 Column Suites Include: Cogent Silica-C™, Bidentate C8™, Bidentate C18™, UDC-Cholesterol™ and Diamond Hydride™		
40100-05P	5 Column Suite	4.6mm x 50mm
40100-05P-2	5 Column Suite	2.1mm x 50mm
40100-05P-3	5 Column Suite	3.0mm x 50mm
40100-7.5P	5 Column Suite	4.6mm x 75mm
40100-7.5P-2	5 Column Suite	2.1mm x 75mm
40100-15P-1	5 Column Suite	1.0mm x 150mm
40100-15P-2	5 Column Suite	2.1mm x 150mm
40100-15P	5 Column Suite	4.6mm x 150mm
40100-25P	5 Column Suite	4.6mm x 250mm

Other sizes and Guard Column suites are available, please inquire. L-76



MicroSolv Technology Corporation
www.MTC-USA.com • 732-578-1777