

Catalog Number

49305-10

Tubing, PEEK, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD,
general grade, solid blue, 5 meter roll, ARE-Applied Research brand. 1
EA.



Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives.

- ID and OD tolerances: +/- 0.05 mm
- Color code: blue
- Pressure rating: 5600 psi (386.1 bar)

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$96.30

Printed On 2024/05/20