

Catalog Number

02001-NAOH

NaOH Conditioning Solution, CE Grade. 0.1N, 100ml. MicroSolv
Brand. 1 EA.



Additional Info

This 0.1 Normal, Sodium Hydroxide Solution is often used for Conditioning CE Capillaries before Injection.

Using premade Conditioning Solutions is very Convenient and can add to Data Precision due to consistency as our SOP.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$37.38

Printed On 2024/05/20