

Catalog Number

46026-FP

Ferrule Puller. Removal Tool for Stainless Steel Ferrules. MicroSolv Brand. 1 EA.



Additional Info

When a ferrule come off stainless steel tubing inside an HPLC column's end fitting or some other port on the instrument, it can be very costly. If one is not careful, you could easily damage the receiving port where the ferrule is stuck.

The Ferrule Puller Tool™ is designed with a "counter" screw thread with sharp edges that will dig into the Ferrule. Then using a counter clock wise turn you can then extract the ferrule without damaging the column or the instrument.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$86.05

Printed On 2024/05/20