

## Catalog Number

49302-CX-09-010

**Union, Connector Kit, Cross**, high pressure, 0.8 mm / 0.031 inch through hole with 1/4-28 screw threads. PEEK, round with an internal "X" flow pattern. Use with 1/16 inch tubing. ARE-Applied Research brand.



### Additional Info

This "Union" is used to connect PEEK or stainless steel tubing the same connectors on all ends in the high pressures zones of HPLC systems. A "Cross (X)" internal flow pattern is often used to split flow into two or three directions or to combine two or three streams into one.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,155.60**

Printed On 2024/05/20