

Catalog Number

9502S-09-PS-M

Caps, closed screw top, solid, with a silicone rubber / polypropylene liner, in "knurled", black caps. AQ Brand. 1000/CS.



Additional Info

Polyfluoroalkyl substances (PFAS) can be tested or stored in our autosampler vials using these caps. The silicone rubber (for optimal sealing), and polypropylene fitted liners (for inertness) are excellent for storage of your sample before or after chromatography. These caps do not contain PTFE. For more information on PFAS, Click [HERE](#).

Use with MicroSolv, AQ™ or RSA™ branded 9mm screw thread autosampler vials for best performance and sealing.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$315.80

Printed On 2024/05/20