

REMOVE KINKS

from Stainless Steel HPLC Tubing for

BETTER FLOW



- Effortlessly remove kinks by placing tubing into the precision hole of the Beta-Tool-2™ HPLC Tubing Straightener and clamping down
- Expertly restore the outer diameter of stainless steel tubing without compromising the inner diameter
- Easily bend tubing without kinking the inside
- Enhance data collection with improved flow

Catalog Number	Product Description
46020-BT	Beta Tool-2. Stainless Steel Tubing straightener with soft grips for tight spaces. Designed for service engineers, used to remove kinks or bends in standard 1/16 In and 1/8 In OD HPLC stainless steel tubing.

MICROSOLV TECHNOLOGY CORPORATION

Leaders in Separation & Purification Technologies

Learn more at www.mtc-usa.com