

# Mobile Phase Bottles

## Mobile Phase Bottles for HPLC

SAFETY COATED

MicroSolv Mobile Phase Borosilicate Glass bottles coated in plastic are a must have for every Laboratory that considers safety their top priority.

The outside of these borosilicate glass bottles are coated in a proprietary plastic that will not fully breach if dropped, even when full with solvent.

- » Safety, Safety, Safety
- » Large Capacities Available
- » Borosilicate Glass Bottles with Plastic Coating on the outside for Safety

They are available with a GL45 neck in various volumes and also, with a Standard Tapered Neck 29/32 in a 1 liter size.

Compatible with most mobile phase filtration devices, you can safely filter directly into these bottles and cap them for storage or transportation.

- » Safety, Safety, Safety
- » Large Capacities Available
- » Borosilicate Glass Bottles with Plastic Coating on the outside for Safety
- » Compatible with Most HPLC Mobile Phase Delivery Systems

## Ordering Information

| Product Number | Volume | Neck Type            | Included Cap Type |
|----------------|--------|----------------------|-------------------|
| 58801-02       | 1L     | Standard Tapered Top | NONE              |
| 58802-01       | 1L     | GL45 Thread          | Solid, Blue       |
| 58819-05       | 2L     | GL45 Thread          | Solid, Blue       |
| 58819-05R      | 2L     | GL45 Thread          | Open Top, Red     |
| 58802-05       | 5L     | GL45 Thread          | Solid, Blue       |
| 58802-10       | 10L    | GL45 Thread          | Solid, Blue       |



## Available Replacement Caps for GL45 Bottles

| Product Number | Delivery Cap Description   | Color |
|----------------|--|-------|
| 58820-06       | Open Hole Screw Cap with a PTFE Insert with 3 Holes, includes PEEK Plug and O-Ring | Blue  |
| 58820-01       | Solid Screw Cap without a liner  | Blue  |
| 58820-10       | Open Hole Screw Cap  | Blue  |

Other Sizes and Formats are Available Upon Request.

MANUFACTURED BY:

**MICROSOLV** TECHNOLOGY CORPORATION

9158 Industrial Blvd NE  
Leland, NC 28451

p: 1.732.380.8900  
f: 1.910.769.9435

customers@mtc-usa.com  
www.mtc-usa.com