

# Crimp Caps

GC or HPLC

## SILICONE RUBBER AT LOW COST • IMPROVE RESEALING • LESS EXTRACTABLES



Standard, aluminum crimp caps for serum or snap-top finish 11mm, 2ml vials. Our caps are made with silicone rubber instead of natural red rubber and therefore offer better resealing, softer on the instrument and needles are cleaner; at the same cost.



Natural red rubber was originally used because it was low cost but now, you can have the advantages of silicone rubber at the cost of or lower than natural rubber.

Available in standard silver, red, green, blue or gold color.

Cat. No.	Description	For more information visit <a href="http://www.mtc-usa.com">www.mtc-usa.com</a>
9502C-40T-B	Cap, crimp, blue aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 100/pack.	
9502C-40T-B-M	Cap, crimp, blue aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 1,000/case.	
9502C-40T-D	Cap, crimp, gold aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 100/pack.	
9502C-40T-D-M	Cap, crimp, gold aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 1,000/case.	
9502C-40T-G	Cap, crimp, green aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 100/pack.	
9502C-40T-G-M	Cap, crimp, green aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 1,000/case.	
9502C-40T-R	Cap, crimp, red aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 100/pack.	
9502C-40T-R-M	Cap, crimp, red aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 1,000/case.	
9502C-40T-S	Cap, crimp, silver aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 100/pack.	
9502C-40T-S-M	Cap, crimp, silver aluminum with fitted red rubber/PTFE septa. Use with snap or crimp top 12 x 32mm autosampler vials. 1,000/case.	

MANUFACTURED BY:

**MICROSOLV** TECHNOLOGY CORPORATION

9158 Industrial Blvd NE  
Leland, NC 28451

p: 1.732.380.8900  
f: 1.910.769.9435

customers@mtc-usa.com  
www.Cogent-HPLC.com