

Vials, Crimp Top - Serum Finish

GC or HPLC

STANDARD 11MM (ND11) SERUM FINISH • EXCELLENT SEALING WITH CRIMP CAPS • WIDE OPENING



Standard crimp top style vials for HPLC or GC autosamplers. The serum finish, 11mm (ND11), allows for crimp caps to form excellent seals. The wide opening of these vials allows for easy filling and decanting and prevents “vapor lock” during filling.

Used in Agilent, CTC, Dani, Fisons, Gerstel, Jasco, Perkin-Elmer, Shimadzu, Spark, Thermo and Varian autosamplers.

Also used for micro-reactions and concentration studies.

Cat. No.

Description

For more information visit www.mtc-usa.com

9502C-0-OCV	Vials, autosampler 2ml 12 x 32mm serum finish/crimp top with wide opening clear borosilicate glass. 100/pack.
9502C-0-WAV	Vials, autosampler 2ml 12 x 32mm serum finish/crimp top with wide opening amber borosilicate glass with write on patch. 100/pack.
9502C-0-WCV	Vials, autosampler 2ml 12 x 32mm serum finish/crimp top with wide opening clear borosilicate glass with write on patch. 100/pack.
9532C-0-WCV	Vials, autosampler 300µl 12 x 32mm serum finish/crimp top with fused insert, clear borosilicate glass with write on patch. 100/pack.

MANUFACTURED BY:

MICROSOLV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

customers@mtc-usa.com
www.Cogent-HPLC.com