

Head Space Vials & Caps

Screw Top for Better Sealing and Convenience



These vials and caps are a convenient and effective alternative to ordinary crimp top, headspace vials and caps. Easy to use, these vials and caps are an extremely reliable system. The caps are more likely to be placed on the vial properly than crimp top vials which require crimping tools.

Built to withstand high internal pressures needed for headspace analysis, these thick wall vials are guaranteed not to burst. Since the caps are the most likely part of the system to leak, these screw caps are made with precision threads for tighter sealing and more reliability.

Autosampler Compatibility:

CTC Pal (Leap)	Shimadzu
Varian	Agilent
Gerstel	Others
Atas	

- **PRECISION THREADS PREVENT LEAKING**
- **AVAILABLE IN DIFFERENT SEPTA TYPES**
- **MAGNETIC CAPS W/O USING INSERTS**
- **HAND TIGHTEN WITH EASE, NO TOOLS REQUIRED**



Cat. No.

Description

For more information visit www.mtc-usa.com

95010-OCV	Vials, Head Space, 10ml, 46 x 22.5mm, screw top, clear borosilicate glass, round bottom. CTC Compliant. 18mm precision threads. 100/pack.
95010-OAV	Vials, Head Space, 10ml, 46 x 22.5mm, screw top, amber borosilicate glass, round bottom. CTC Compliant. 18mm precision threads. 100/pack.
95020-OCV	Vials, Head Space, 20ml, 75.5 x 22.5mm, screw top, clear borosilicate glass, round bottom. 18mm precision threads. 100/pack.
95020-OAV	Vials, Head Space, 20ml, 75.5 x 22.5mm, screw top, amber borosilicate glass, round bottom. 18mm precision threads. 100/pack.
95020-40T-S	Caps, Head Space, silver screw caps, magnetic, 8mm center hole, butyl red/PTFE grey, 55 Shore A, 1.6mm thick septa. 18mm precision threads. 100/pack.
95020-45S-T	Caps, Head Space, silver screw caps, magnetic, 8mm center hole, silicone transparent blue/PTFE white, 45 Shore A, 1.3mm thick septa. 18mm precision threads. 100/pack.
95020-45S-W	Caps, Head Space, silver screw caps, magnetic, 8mm center hole, silicone white/PTFE red, 45 Shore A, 1.3mm thick septa. 18mm precision threads. 100/pack.
95020-55S-W	Caps, Head Space, silver screw caps, magnetic, 8mm center hole, silicone white/PTFE blue, 55 Shore A, 1.5mm thick septa. Temperature resistant from -60°C up to 200°C. 18mm precision threads. 100/pack.
95020-50S-W	Caps, Head Space, silver screw caps, magnetic, 8mm center hole, silicone white/PTFE blue, 50 Shore A, 1.5mm thick septa. 18mm precision threads. 100/pack.

MANUFACTURED BY:

MICRO SOLV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

customers@mtc-usa.com
www.Cogent-HPLC.com