

MRQ™ Sample Vial

More than an AutoSampler Vial



This autosampler vial is completely compatible with most current commercially available autosamplers. Purchase vials individually or in cost saving Easy Purchase packs with caps.

These vials are very useful when you want to recover all of your sample or have only limited volume to start with.

The all glass patent pending construction is produced without any internal “ridges or shelves” for maximum drainage to bottom well.

Unlike vials with fused inserts or inserts with attached springs, these rugged limited volume vials can be centrifuged, lyophilized or used with “finger vortexers”. MRQ vials have a precision formed conical center well in an all borosilicate, clear glass vial.

Some Useful Applications:

- Autosampler Vials
- Library Storage Vials
- Lyophilization Containers
- Centrifugation of Very Small Volumes

- AVAILABLE AS 9 MM SCREW TOP OR SNAP/CRIMP TOP
- 2µL RESIDUAL VOLUME WITH 1.2 ML TOTAL VOLUME
- COMPARE TO “TOTAL RECOVERY VIALS™”
- MAXIMUM RECOVERY FOR SIDE PORT NEEDLES

Cat. No.

Description

For more information visit www.mtc-usa.com

9512S-OCV-T	MRQ™ Max Recovery Vials 1.2ml, 9mm, 12 x 32mm, Screw Top. Caps sold separately. 100/pack.
9512S-OCV-T-RS	MRQ™ Max Recovery Vials 1.2ml, 9mm, 12 x 32mm, Screw Top. RSA™ Glass. Caps sold separately 100/pack.
9512C-OCV-T	MRQ™ Max Recovery Vials 1.2ml, 12 x 32mm, Snap Top. Caps sold separately. 100/pack.
9512C-OCV-T-RS	MRQ™ Max Recovery Vials 1.2ml, 12 x 32mm, Snap Top RSA™ Glass. Caps sold separately 100/pack.

MANUFACTURED BY:

MICROSOLV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

customers@mtc-usa.com
www.Cogent-HPLC.com