Autosampler Vial Racks

For 12 x 32mm Crimp, Snap and Screw Neck Vials





These Autosampler Vial Racks are a safe and convenient way to store and transport your sample vials. The cavities are 130mm x 130mm x 45mm, designed to fit 12mm x 32mm Crimp, Snap and Screw Neck Vials. They are available in 6 assorted colors: blue, orange, yellow, pink, white and green. The cavities allow for your vials to stand safely in a standing position on the laboratory table. They are coded alpha numerically for clear sample identification.

These racks are ideal for space saving storage in cooling units. As the transparent lid prevents condensation on the closures thus avoiding a possible contamination. They also are temperature resistant from -80°C to +100°C.

- CONVENIENT AND SAFE
- FOR STORAGE AND TRANSPORTATION
- STACKABLE
- AUTOCLAVABLE

Cat. No.	Description	For more information visit www.mtc-usa.com
9502SB-01	Storage Box/Rack with cover, PP, for 12 x 32mm autosampler vials, blue, 130mm x 130mm x 45mm. 81 cavities with alpha numeric coding of all 4 sides as well as the cavity bottom. 1 each.	
9502SB-02	- · · · · · · · · · · · · · · · · · · ·	P, for 12 x 32mm autosampler vials, orange, 130mm th alpha numeric coding of all 4 sides as well as the
9502SB-03		P, for 12 x 32mm autosampler vials, pink, 130mm x alpha numeric coding of all 4 sides as well as the
9502SB-04	- · · · · · · · · · · · · · · · · · · ·	P, for 12 x 32mm autosampler vials, yellow, 130mm th alpha numeric coding of all 4 sides as well as the
9502SB-05		P, for 12 x 32mm autosampler vials, green, 130mm x alpha numeric coding of all 4 sides as well as the
9502SB-06	- · · · · · · · · · · · · · · · · · · ·	P, for 12 x 32mm autosampler vials, transparent, vities with alpha numeric coding of all 4 sides as well