


Autosampler Vial Racks

For 15 x 45mm, 4ml Vials

CONVENIENT AND SAFE • FOR STORAGE AND TRANSPORTATION
STACKABLE • AUTOCLAVABLE



These Autosampler Vial Racks are a safe and convenient way to store and transport your sample vials. Temperature resistant from -80°C to +100°C, these special racks are ideal for freezer storage. The transparent lid prevents condensation on the closures to minimize or prevent problems that occur with preparing samples the day before analysis.

The cavities in these racks are 130mm x 130mm x 52mm and designed to securely fit 15mm x 45mm snap and screw top vials. Available in red.



Also available for Headspace Vials
Please inquire

Cat. No.

9504SB-03

Description

Storage Box/Rack with cover, PP, for 15 x 45mm autosampler vials, red, 130mm x 130mm x 52mm. 49 cavities with alpha numeric coding on the cavity bottom. 1 each.

For more information visit www.mtc-usa.com

MANUFACTURED BY:

MICROSOLV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE
Leland, NC 28451

p: 1.732.380.8900
f: 1.910.769.9435

customers@mtc-usa.com
www.Cogent-HPLC.com