



Metal Content Affects Silica Particle and Contributes to Peak Shape - Tech Information

Date: 22-APRIL-2012 Last Updated: 1-DECEMBER-2025

The electron withdrawing properties of trace metals in standard HPLC silica particles enhances the negativity of surface silanol groups and contributes to peak tailing when high levels of trace metals are present.

NOTE: This much less of a concern for TYPE-C™ silica.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com