



Precision and Suitability of the Cogent TYPE-C HPLC Columns in LCMS - Tech Information

Date: 14-APRIL-2020 Last Updated: 30-NOVEMBER-2025

The Cogent TYPE-C™ Silica HPLC columns produce lower background noise in LCMS than most other HPLC columns. For column precision, typically %RSD for retention time of replicate injections is less than 1%.

As for detection precision, we highly recommend using an appropriate internal standard in LCMS. This will help with reliable quantitation of your analytes.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com