

## What should I know about Sample Loops - FAQ

*Date: 24-MAY-2014 Last Updated: 18-JUNE-2025*

Sample loops are available in a range of internal volumes to suit different injection needs in HPLC and other analytical systems.

Depending on the volume:

- Some loops will have multiple visible turns of tubing.
- Others, especially those with very small volumes, may appear straight or without a visible loop.

This is completely normal. In cases where the tubing is too short to form a physical loop, the product still functions exactly as intended. The presence or absence of a visible loop does not affect performance or accuracy.

---

 Click [HERE](#) for sample loop ordering information and product images



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)