

## Understanding 316 Stainless Steel Composition and Applications - HPLC Primer

*Date: 12-MAY-2012 Last Updated: 19-JUNE-2025*

Grade 316 stainless steel is a high-quality, corrosion-resistant alloy widely used across industries, including laboratory environments. At MicroSolv, we utilize 316 stainless steel in our laboratory tubing, fittings, and HPLC columns due to its superior performance characteristics.

### What Is Grade 316 Stainless Steel?

According to Atlas Steels Australia, Grade 316 is the standard molybdenum-bearing grade among austenitic stainless steels. It is second only to Grade 304 in commercial importance. The addition of molybdenum enhances its resistance to:

- Pitting and crevice corrosion, especially in chloride-rich environments
- General corrosion, making it ideal for harsh chemical and marine conditions

### Key Properties of 316 Stainless Steel:

- Excellent corrosion resistance in a wide range of environments
- Superior formability and weldability, suitable for complex shapes and assemblies
- No post-weld annealing required for thin sections
- High strength and durability, even at elevated temperatures

These properties make Grade 316 stainless steel a preferred material in industrial, architectural, transportation, and laboratory applications.

### Additional Notes on Material Designations:

- Tempering is a heat treatment process used to increase the toughness of iron-based alloys. It differs from titanium designations such as Ti1 or Ti2, which refer to the percentage of pure titanium in an alloy. Generally, a higher number indicates a greater titanium content and improved performance in specific applications.

---

 Click [HERE](#) to explore our Stainless Steel Laboratory Tubing