

## Why is it Not Recommended to Use FEP Tubing Above 80°C When the Melting Point is Much Higher - FAQ

*Date: 26-OCTOBER-2015 Last Updated: 8-NOVEMBER-2025*

**FEP polymer has a melting point of 260°C and a general recommended maximum operating temperature during usage of around 200°C.**

If pressure is applied in the tubing, however, the maximum temperature should be lowered accordingly. The maximum temperature during usage should be 80°C under these conditions. If no pressure is applied, the tubing can be used up to 200°C.

Click [HERE](#) for low pressure HPLC Tubing ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)