



Standard silicone pre-slit sealing mat versus solvent resistant mat - INTERNAL ONLY

Silicone rubber mats are only suitable with pure aqueous and buffer solutions.

The solvent resistant mat should be used with Methanol, Acetonitrile and water/buffers. These should not be used with non-polar solvents or aggressive acid or base.

U-2D Silicone Mats Product Page

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com