

Do You Have Autosampler Vial Screw Caps Without Septa - FAQ

Yes! We offer two types of caps without septa:

The first type is a completely closed cap where the area that would normally have a hole for the septum is covered. These are offered with **liners** instead of septa. Please see catalog numbers [9502S-09-SR](#) and [9502S-09-SR-K](#) for examples.



The second type is simply the polypropylene cap in which the septa is **Not present**. The caps are normally referred to as open hole caps. An example catalog number for this type of cap is [9502S-00-B-M](#).

