

Which Grade of Stainless Steel Tubing and Column Hardware do you Sell - FAQ

Date: 2-JULY-2013 Last Updated: 26-OCTOBER-2025

The alloy used in our stainless steel tubing and our HPLC columns is SS316.

This is also true for filters/sparging stones.

Click [here](#) to learn more about what this nomenclature means in terms of stainless steel grade.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com