



Are the Last Drop Inlet Filters UHPLC Compatible - FAQ

Date: 22-JUNE-2020 Last Updated: 8-NOVEMBER-2025

YES.

MICROSOLV's Last Drop™ inlet filters for mobile phase bottles ARE compatible with UHPLC, but it should be noted that for columns with 3µm or smaller particle size packing including UHPLC columns, a 0.2µm pore size filter is recommended. *Flow will be slower due to the smaller pores in the filter.*

Click [HERE](#) for Last Drop Filters™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com