

What is Best Way to Know if There is a Crack in a CE Capillary - FAQ

Date: 28-APRIL-2012 Last Updated: 3-NOVEMBER-2025

The best way to know is to look at the entire capillary under either a stereo microscope or using a hand held magnifier offered on our website.

Click [HERE](#) for ordering information and pictures of the magnifier for this purpose.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com