

What Are the Dimensions of 4.6mm ID 5um Frits used with matching PEEK Column Hardware - FAQ

Date: 25-NOVEMBER-2020 Last Updated: 28-NOVEMBER-2025

The Frit itself has an outer diameter of 0.188" and a thickness of 0.062" and is in a PEEK ring with 0.25" O.D. and 0.181" ID. The Frit is press-fit into the ring without glues or binders.

Click [HERE](#) for PEEK column hardware ordering page.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com