

Keep Surface Treated Frits From Biofouling - Tips & Suggestions

Date: 11-SEPTEMBER-2021 Last Updated: 28-NOVEMBER-2025

Use non surfactant, Brij-35 to reduce biofouling.

When many proteins get trapped on or adhere to stainless steel frits that are hydrophobic like our metal free surfaces, it is a good idea to add Brig-35 to a wash buffer to facilitate removal of very sticky proteins including BSA, mouse IgG and human plasma proteins.

Brig-35 will not diminish or change the effectiveness of the metal free surface. It is important to work in low concentrations of any ionic polymer to minimize biofouling due to surfaces.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com