



## What is the best way to track samples in U-2D inserts - Tips and Suggestions

The U-2D™ inserts can be written on using our laboratory markers found [HERE](#).

Also, the U-2D™ racks have positions that are marked with alpha-numeric embossed markings on the sides of the racks.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)