


## Fluvoxamine Maleate USP Assay with HPLC – AppNote

### Alternate Column Used for a Similar USP Assay Method of Fluvoxamine Maleate


This similar USP fluvoxamine maleate Assay is easily performed, and the 5 injections overlay below demonstrates run to run consistency and separation of E/Z isomers. The peak efficiency is very good and peak tailing guidelines are easily met with RSD values of less than 1.0%.

This demonstrates a good alternate column for your fluvoxamine maleate USP method.



#### Peak:

Fluvoxamine maleate (E/Z isomers)



#### Met hod Conditions:

**Column:** Cogent RP C8™, 5um, 100Å

**Catalog No.:** [68508-15P](#)

**Dimensions:** 4.6mm x 150mm

---

**Note 1:** For the full method conditions and details, consult the official United States Pharmacopeia–National Formulary (USP–NF.)

**Note 2:** Fluvoxamine maleate tablets are indicated for the treatment of obsessions and compulsions in patients with obsessive compulsive disorder (OCD).

---



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)