

What is the bending radius of 0.0025 inch ID 1/16 inch OD PEEK tubing - FAQ

The smallest bend for PEEK tubing is a radius of 25 mm which is a diameter of 5 cm (2 inch).

Click [HERE](#) for the ordering information page for PEEK tubing.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com