



## FlexiTube specifications - Tech Information

### The FlexiTube™ length and coil data for selecting the right tube for you.

Ordering information for each type of FlexiTube can be found on the [ordering information](#) page, under each individual item.

| Type | Coil resting position | Coil extended position | Coil quantity | FlexiTube length fully extended | Total length of tubing |
|------|-----------------------|------------------------|---------------|---------------------------------|------------------------|
| 2    | 10 mm                 | 30 mm                  | 4             | 330 mm                          | 710 mm                 |
| 3    | 20 mm                 | 120 mm                 | 9             | 420 mm                          | 945 mm                 |
| 4    | 35 mm                 | 200 mm                 | 15            | 500 mm                          | 1225 mm                |
| 5    | 50 mm                 | 300 mm                 | 23            | 600 mm                          | 1585 mm                |
| 6    | 100 mm                | 400 mm                 | 46            | 700 mm                          | 2060 mm                |
| 7    | 150 mm                | 500 mm                 | 69            | 900 mm                          | 3150 mm                |
| 8    | 170 mm                | 8000 mm                | 78            | 1100 mm                         | 4710 mm                |
| 9    | 225 mm                | 1100 mm                | 104           | 1400 mm                         | 6100 mm                |
| 10   | 250 mm                | 1300 mm                | 116           | 1600 mm                         | 6930 mm                |

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)