



Are the MICROSOLV Mobile Phase Bottles Type 1, Class A Glass - FAQ

Date: 9-FEBRUARY-2023 Last Updated: 8-NOVEMBER-2025

Yes, they are.

Both the safety-coated and non-coated MICROSOLV mobile phase bottles are crafted from high-quality borosilicate 3.3 glass, a type of Type 1, Class A borosilicate glass. This material offers excellent chemical resistance, making these bottles suitable for a wide range of laboratory applications.

- Click [HERE](#) for ordering information for safety coated bottles.
- Click [HERE](#) for ordering information for non - coated bottles.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com