



Physical Nature of PEEK Frits Caution - Technical Information

Date: 6-NOVEMBER-2023 Last Updated: 28-NOVEMBER-2025

Due to the physical nature of the material and how it is manufactured, one must be careful when handling PEEK Frits.

PEEK frits can be brittle. Do not throw the bag of frits or drop it. Use caution when removing the frits in the bag, as bumping into one another can cause dislodging of the material.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com